

ATOMO Dental, Inc.

Add: 1241 Quarry Lane #125, Pleasanton, CA 94566

Tel: 866-9800988 / 510-3642700

Email: support@atomodental.com

Web: www.atomodental.com

Blue Nitrile Examination Gloves Product Specification

Type	Nitrile Exam Glove, Powder-Free, Non-Sterile, 3.5g													
Grade	Medical			Surface	Finger textured				Packing	100pcs/box,10boxes/case (packaging can be customized)				
Product	Color	Size	Weight (g)	Physical Dimension					Physical Property			Test AQL	Certification	
				Length (mm)	Palm Width (mm)	Finger Thickness (mm)	Palm Thickness (mm)	Cuff Thickness (mm)	Tensile Strength (MPa)	Elongation (%)				
										Before Aging	After Aging			
Powder Free Nitrile Exam Gloves	Blue	XS	2.5±0.2	240±5	75±5	>0.08	>0.07	>0.05	>14	500	400	AQL1.5	510K, CE certification, ISO9001, ISO13485	
		S	3.0±0.2		85±5									
		M	3.5±0.2		95±5									
		L	4.0±0.2		105±5									
		XL	4.5±0.2		115±5									
Our manufacturer has established the Quality management system in according with ISO9001 and ISO 13485 standards														
Sampling Procedures according to the standard of ISO2859, ASTM D6319, Inspection Level for physical dimensions and force / elongation at break is according to S-2, water test according to G-1														